




# ABOUT ISOLAB


ISOLAB is a high quality laboratory consumables manufacturer, especially well-known for its excellent glassware products, and today preferred by thousands of lab professionals in more than 110 countries all around the World.

ISOLAB Laborgeräte GmbH is founded in 1997, in Wertheim Germany. Today ISOLAB GmbH's main office and Europe Warehouse is located in Eschau, Bayern.

 Main Office & Warehouse  
Eschau, Bayern - Germany




Exported to  
110 COUNTRIES

 11.190 m<sup>2</sup> space  
WAREHOUSE



266 DISTRIBUTORS  
in 110 COUNTRIES

 A product portfolio of  
3.584 PRODUCTS



2 EXPORT BASES  
Eschau – European Sales  
Istanbul – Rest of the World



Qualified with  
DIN EN ISO 9001



ALL  
FOR  
LAB

ISOLAB

glasswares   
consumables   
equipments   
instruments   
chemicals



## ABOUT ISOLAB

Starting on its foundation in 1997, ISOLAB GmbH has been a company that relies on its products of highest quality, puts great emphasis on customer satisfaction as well as service quality and sustains the strength of its past into the future.

Based on the principles of innovative approach and entrepreneurial vision, ISOLAB has been dedicated to provide latest products and solutions to its customers for daily work in universities, institutions, laboratories, research and applications in various industries. The success we achieved in laboratory consumables is a clear evidence of that customer trust.

By carrying out the correct moves in accordance with the strategies we have developed, we as the ISOLAB Company, have transformed into a global company and achieved a great deal of success, both at home and internationally. During this period, when we are daily consolidating our leadership and trying to keep abreast of the pace with which the world is changing, we are inevitably focusing on investing time and resources in our brand.

Combining the strengths of premium quality with a visionary spirit, enriching these with extensive know-how of almost 25 years and also offering all of this at a fair price, ISOLAB positions itself as a leading company and a very well-known brand in the world.



## PRODUCT LINE

- Volumetric Line & Liquid Handling
- LAB CONSUMABLES & LABWARE
- Distillation, Filtration Separation & Sample Preparation
- Cleaning & Sterilization Protection & Safety
- Density, Sedimentation, Temperature, Time pH & Weight Measurement
- Stirring, Mixing, Shaking, Heating & Cooling
- CHROMATOGRAPHY & SPECTROSCOPY
- Life-Science



## CERTIFICATIONS

ISOLAB products are being produced in accordance with many international quality standards.

→ISO 17025

→ISO 13485

→ISO 9001

→ISO 14001

→ISO 45001

→FOOD CONTACT CONFORMITY

## TOP GLASS PRODUCTS

→ Glass flasks are manufactured from heat resistant borosilicate 3.3.


→ Borosilicate 3.3 glass flasks are manufactured according to DIN EN ISO 1042 standard.

→ Calibrated according to DIN/ISO Class "A" standard.

Manufactured from heat resistant borosilicate 3.3 glass according to DIN EN ISO 4788 standard. →

Calibrated by automatic robots according to DIN/ISO Class "A" standard. →


Conformity mark ("DE-M") and batch number printed on the body surface. →



→ Manufactured from heat resistant, durable borosilicate glass according to DIN 12576 and DIN 12581 standards. →

→ Featured with glass side arms fitting 8-9 mm bore size tubings. →

→ Thick wall, and reinforced rim prevent breakage during washing and transportation. →







## SOME OF PLASTIC PRODUCTS



→ Manufactured from ultra clear polypropylene according to ISO 7056 standard with thick walls.

→ Resistant to most acids, bases and many commonly used solvents.

→ Functional large spout minimize possible spills during liquid transfer.

Manufactured from ultra clear, high grade polypropylene. →

Ultra rugged tube wall guarantees resistance to high centrifugal forces up to 6.000 RCF. →

Autoclavable at 121°C for 20 minutes. →

Caps should not be kept closed during autoclaving process. →

